


Using a microscope

How do you use a microscope to help see onion cells clearly?



Statement list

<p>Look down the eyepiece lens.</p>	<p>Turn the knob upwards to bring the onion cells into focus.</p>
<p>Put them on a clean microscope slide with a drop of water.</p>	<p>Put the slide onto the stage of the microscope.</p>
<p>Turn the knob so the second lens almost touches the slide.</p>	<p>Peel a thin layer of cells from the inner part of an onion.</p>
<p>Put them in water to stop them from drying out.</p>	<p>Switch on the light to illuminate the onion cells.</p>